NADAL® COVID-19 S1-NAb Test (test cassette)

REF 243301N-20

de	Gebrauchsanweisung	2	pt	Instruções de Utilizaçāo	26
en	Instructions for use	6	cs	Návod k použití	30
fr	Instructions d'utilisation	10	fi	Käyttöohje	34
es	Instrucciones de uso	14	sv	Användarinstruktioner	38
it	Istruzioni per l'uso	18		Symbols	43
pl	Sposób użycia	22		Our Teams	44

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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15. Referenzen

- Coronavirus disease 2019 (COVID-19) in the EU/EEA and the UK – ninth update, 23 April 2020. Stockholm: ECDC; 2020.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81:85-164.
- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. TrendsMicrobiol 2016;24:490-502.

1. Intended Use

The NADAL® COVID-19 S1-NAb Test is a lateral flow chromatographic immunoassay for the qualitative detection of potentially neutralising antibodies to SARS-CoV-2 in human whole blood, serum or plasma specimens (see section 12 'Limitations'). This test is intended for use as an aid in determining whether an immune response to the S1 protein of SARS-CoV-2 has occurred. The test procedure is not automated and requires no special training or qualification. The NADAL® COVID-19 S1-NAb Test is designed for professional use only.

2. Introduction and Clinical Significance

COVID-19 (Corona Virus Disease) is an infectious disease caused by the recently discovered coronavirus SARS-CoV-2. The most common symptoms of COVID-19 are fever, a dry cough, fatigue, sputum production, shortness of breath, a sore throat and headaches. Some patients may have myalgia, chills, nausea, nasal congestion and diarrhoea. These symptoms begin gradually and are mild in most cases. Some people become infected but do not develop any symptoms nor feel unwell. Most people (about 80%) recover from the disease without special treatment. Approximately one in six people who get infected with COVID-19 become seriously ill and develop difficulty breathing. Elderly people, and those with pre-existing conditions, such as high blood pressure, heart problems or diabetes, are more likely to develop serious illness.

COVID-19 is transmitted via respiratory droplets that are exhaled by infected people via coughing, sneezing or talking. These droplets can be inhaled or ingested directly by other people or can contaminate surfaces, which can then remain infectious for several days. Most estimates of the incubation period for COVID-19 range from 1 to 14 days, during which time people might already be infectious without displaying disease symptoms.

3. Test Principle

The NADAL® COVID-19 S1-NAb Test is a lateral flow chromatographic immunoassay for the qualitative detection of potentially neutralising antibodies to SARS-CoV-2 in human whole blood, serum or plasma specimens.

Anti-IgM and anti-IgG antibodies are pre-coated onto the test line region (T) of the membrane. During testing, the specimen reacts with spike protein antigens which are conjugated to coloured particles. The mixture then migrates along the membrane chromatographically by capillary action and reacts with anti-IgM and anti-IgG antibodies in the test line region (T) of the membrane. The presence of a red line in the test line region (T) indicates a positive result. The absence of a red line in the test line region (T) indicates a negative result.

The colour change from a blue line to a red line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® COVID-19 S1-NAb test cassettes*
- 20 capillary tubes
- 1 buffer (3 mL)**

- 1 package insert

*containing the preservative sodium azide: <0.02% (7.5 ng/test)

**Phosphate buffer containing the following preservatives: sodium azide: 0.2 mg/mL and kanamycin sulfate: 0.25 g/L

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP, Concentrations are below the exemption threshold.

5. Additional Materials Required

- Specimen collection containers (suitable for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible

pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).

- Temperature can adversely affect test results.
- Failure to bring specimens and reagents to room temperature prior to testing may decrease assay sensitivity. Inaccurate or inappropriate specimen collection, storage or transport may yield false negative test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

The NADAL® COVID-19 S1-NAb Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as K₂EDTA, sodium citrate, sodium heparin or lithium heparin should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 2 days of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.


If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
 2. Place the test cassette on a clean and level surface.
 3. Holding the capillary tube vertically, draw the specimen (whole blood/serum/plasma) up to the fill line (approximately 10 µL) and add it to the specimen well (S) of the test cassette.
- Alternatively, a micropipette (10 µL) may be used.**
4. Holding the buffer bottle vertically, add 2 drops of buffer (approximately 80 µL) to the buffer well (B). **Avoid air bubbles forming.**
 5. Start the timer.
 6. Wait for the red line(s) to appear. Read the test result after exactly 10 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive

The blue line in the control line region (C) turns red. Another red line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the concentration of neutralising antibodies to SARS-CoV-2 in the specimen. Therefore, any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


Negative

The blue line in the control line region (C) turns red. No line develops in the test line region (T).

Invalid

Results from any test where the blue line in the control line region (C) remains completely or partially blue at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette: The line turning from blue to red in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® COVID-19 S1-NAb Test is for professional *in-vitro* diagnostic use only. It should be used for the qualitative detection of potentially neutralising antibodies to SARS-CoV-2 in human whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in the concentration of potentially neutralising antibodies to SARS-CoV-2 can be determined with this qualitative test.
- The NADAL® COVID-19 S1-NAb Test only detects the presence of potentially neutralising antibodies to SARS-CoV-2 in specimens and should not be used as the sole criterion for determining whether an immune response to the S1 protein of SARS-CoV-2 has occurred.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- The continued presence or absence of antibodies cannot be used to determine the success or failure of treatment or vaccination.
- Results from immunosuppressed patients should be interpreted with caution.
- A high-dose hook effect may occur where the colour intensity of test lines decreases as the concentration of neutralising antibodies to SARS-CoV-2 increases. If a high-dose hook effect is suspected, the dilution of specimens may increase the colour intensity of test lines.

13. Expected Values

The response of the human immune system to the S1 protein of SARS-CoV-2 is connected to the formation of specific antibodies against different epitopes of this protein sub-unit. Only some of the IgM and IgG antibodies formed have a neutralising effect on SARS-CoV-2 virions. IgG antibodies can indicate immunity. In some cases, it can be that antibodies are first present in detectable numbers quite some time after the infection or vaccination has occurred. In cases of a sufficient innate or a sufficient cellular immune response, an absence of detectable antibodies is also possible.

14. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® COVID-19 S1-NAb Test was evaluated in comparison with a microneutralisation assay (cut-off titer 1:20) using clinical specimens from recovered patients as well as healthy, unvaccinated individuals.

NADAL® COVID-19 S1-NAb Test	Microneutralisation assay		
	Positive	Negative	Total
	30	0	30
	0	18	18
Total	30	18	48

Diagnostic sensitivity: >99.9% (90.5% - 100.0%)*

Diagnostic specificity: >99.9% (84.7% - 100.0%)*

Overall agreement: >99.9% (93.9% - 100.0%)*

*95% confidence interval

The NADAL® COVID-19 S1-NAb Test was evaluated in comparison with a microneutralisation assay (cut-off titer 1:20) using clinical specimens from vaccinated as well as healthy, unvaccinated individuals.

NADAL® COVID-19 S1-NAb Test	Microneutralisation assay		
	Positive	Negative	Total
	30	2**	32
	0	18	18
Total	30	20	50

Diagnostic sensitivity: >99.9% (91.1% - 100.0%)*

Diagnostic specificity: 90.0% (68.3% - 98.8%)*

Overall agreement: 96.0% (86.3% - 99.5%)*

*95% confidence interval

**Both false positive specimens contained low antibody concentrations and were identified as positive using the microneutralisation assay at a cut-off titer of 1:10.

In a seroconversion study, 15 test subjects were tested using the NADAL® COVID-19 S1-NAb Test before vaccination, 28 days after the first vaccination and 12 days after the second vaccination with an inactivated SARS-CoV-2 vaccine. Before the first vaccination, all the test subjects were negative. Potentially neutralising antibodies to the S1 protein of SARS-CoV-2 could be detected using the NADAL® COVID-19 S1-NAb Test in 5 test subjects after the first vaccination and in all 15 test subjects after the second vaccination. All specimens tested were confirmed using the microneutralisation assay. The specificity of the NADAL® COVID-19 S1-NAb Test was evaluated with 208 serum specimens from June 2019 (prior to COVID-19 outbreak). The NADAL® COVID-19 S1-NAb Test identified 207 out of 208 specimens as true negative and 1 specimen as false positive (specificity 99.5%).

NADAL® COVID-19 S1-NAb Test	Serum specimens from June 2019		
	Positive	1	
	Negative	207	
	Total	208	

Specificity: 99.5% (97.5% - 99.9%)*

*95% confidence interval

The NADAL® COVID-19 S1-NAb Test was evaluated in comparison with an ELISA method using clinical specimens from vaccinated as well as healthy, unvaccinated individuals.

NADAL® COVID-19 S1-NAb Test	ELISA		
	Positive	Negative	Total
	30	0	30
	0	30	30
Total	30	30	60

Diagnostic sensitivity: >99.9% (90.5% - 100.0%)*

Diagnostic specificity: >99.9% (90.5% - 100.0%)*

Overall agreement: >99.9% (95.1% - 100.0%)*

*95% confidence interval

Analytical performance

Detection limit

The detection limit of the NADAL® COVID-19 S1-NAb Test is 50 ng/mL for neutralising antibodies to SARS-CoV-2.

Rev. 1, 2021-03-02 OM

Measuring range

No adverse effect on T-line formation (prozone effect) was observed when testing specimens containing a concentration of neutralising antibodies to SARS-CoV-2 as high as 1 mg/mL. Thus, the measuring range of the test is between at least 50 ng/mL and 1 mg/mL.

Interfering substances

Solutions of the following potentially interfering substances were tested in 2 test series (without neutralising antibodies to SARS-CoV-2 and spiked with neutralising antibodies to SARS-CoV-2) and showed no interference with the NADAL® COVID-19 S1-NAb Test at the concentrations listed below.

Acetylsalicylic acid	200 mg/L
Albumin	20 g/L
Acetaminophen	200 mg/L
Ascorbic acid	20 g/L
Bilirubin	10 g/L
Caffeine	200 mg/L
Ethanol	1%
Creatine	2 g/L
Gentistic acid	200 mg/L
Haemoglobin	10 g/L
Oxalic acid	600 mg/L
Uric acid	20 g/L
Triglycerides	600 mg/L

Cross-reactivity

Anti-Influenza A virus, anti-Influenza B virus, anti-RSV, anti-adenovirus, HBsAb, anti-syphilis, anti-*H. pylori*, anti-HIV, anti-HCV, HAMA, MERS-CoV and rheumatoid factors positive specimens were tested using the NADAL® COVID-19 S1-NAb Test. No cross-reactivity with the specimens was observed when tested using the NADAL® COVID-19 S1-NAb Test. Anti-SARS-CoV positive specimens show cross-reactivity with the NADAL® COVID-19 S1-NAb Test.

Precision

Repeatability

Precision was established by testing 10 replicates of negative specimens and specimens positive for neutralising antibodies to SARS-CoV-2.

Reproducibility

Reproducibility was established by testing 5 replicates of a positive control, 8 negative serum specimens and 5 negative whole blood specimens. Testing was performed by 3 operators using 3 independent NADAL® COVID-19 S1-NAb test lots at 3 different sites on 3 separate days.

The NADAL® COVID-19 S1-NAb Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

15. References

- Coronavirus disease 2019 (COVID-19) in the EU/EEA and the UK – ninth update, 23 April 2020. Stockholm: ECDC; 2020.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res* 2011; 81:85-164.
- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol* 2016;24:490-502.

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nal von minden GmbH