# **NADAL® COVID-19 Ag+Influenza A/B Test (test cassette)**

**REF 243203N-20**

<b>DE</b>	Gebrauchsanweisung	2	<b>PT</b>	Instruções de Utilizaçāo	38
<b>EN</b>	Instructions for use	8	<b>CZ</b>	Návod k použití	44
<b>FR</b>	Instructions d'utilisation	14	<b>FI</b>	Käyttöohje	50
<b>ES</b>	Instrucciones de uso	20	<b>SE</b>	Användarinstruktioner	56
<b>IT</b>	Istruzioni per l'uso	26		Symbols	63
<b>PL</b>	Sposób użycia	32		Our Teams	64

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Carl-Zeiss-Strasse 12  
47445 Moers  
Germany

Moers  
Tel: +49 (2841) 99820-0  
Fax: +49 (2841) 99820-1

Regensburg  
Tel: +49 941 29010-0  
Fax: +49 941 29010-50

[www.nal-vonminden.com](http://www.nal-vonminden.com)  
[info@nal-vonminden.com](mailto:info@nal-vonminden.com)

Directors:  
Sandra von Minden  
Roland Meißner  
Thomas Zander

Commercial reg. Kleve  
HRB 5679  
Steuer-Nr. 244/133/00130  
UST-ID-Nr. DE 189 016 086

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### 15. Referenzen

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

Rev. 1, 2020-10-12 OM

## 1. Intended Use

The NADAL® COVID-19 Ag+Influenza A/B Test is a lateral flow chromatographic immunoassay for the simultaneous qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens (see section 12 'Limitations'). This test is intended for use as an aid in the diagnosis of infections with SARS-CoV-2 as well as influenza virus type A and B. Note that the concentration of viral nucleoprotein antigens may vary in the course of the disease and might fall below the detection limit of the test. Possible infectiousness of test subjects cannot be ruled out based on negative test results. The test procedure is not automated and requires no special training or qualification. The NADAL® COVID-19 Ag+Influenza A/B Test is designed for professional use only.

## 2. Introduction and Clinical Significance

COVID-19 (Corona Virus Disease) is the infectious disease caused by the recently discovered coronavirus SARS-CoV-2. The most common symptoms of COVID-19 are fever, dry cough, fatigue, sputum production, shortness of breath, sore throat and headache. Some patients may have myalgia, chills, nausea, nasal congestion and diarrhoea. These symptoms begin gradually and are mild in most of the cases. Some people become infected but do not develop any symptoms and do not feel unwell. Most people (about 80%) recover from the disease without special treatment. Approximately one in six people who get infected with COVID-19 becomes seriously ill and develops difficulty breathing. Elderly people, and those with pre-existing conditions, such as high blood pressure, heart problems or diabetes, are more likely to develop serious illness. COVID-19 is transmitted via respiratory droplets that are exhaled by infected people via coughing, sneezing or talking. Most estimates of the incubation period for COVID-19 range from 1 to 14 days, during which people might already be infectious without showing disease symptoms.

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, influenza A subtypes are differentiated by the antigen variability of surface glycoproteins (haemagglutinin and neuraminidase). The influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, the elderly or those with chronic respiratory diseases. However, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections.

With a sensitivity and specificity of almost 100%, RT-PCR is the gold standard of laboratory diagnostics. However, as a laboratory method, it does not allow immediate diagnosis and requires qualified personnel and expensive laboratory equipment. Antiviral influenza drugs, which have been available since the late 1990s, are most effective if administered soon after the onset of the illness (within 48 hours). Rapid tests for the detection of influenza antigens can

therefore enable successful antiviral therapy to be carried out in good time, thereby improving patient healthcare and reducing health costs.

## 3. Test Principle

The NADAL® COVID-19 Ag+Influenza A/B Test is a lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens.

The NADAL® COVID-19 Ag+Influenza A/B Test contains two internal test strips (left side of the test cassette: COVID-19 Ag Test; right side of the test cassette: Influenza A/B Test).

Anti-SARS-CoV-2 antibodies are immobilised in the test line region 'T' of the membrane. Anti-influenza virus type A and B antibodies are immobilised in the test line regions 'A' and 'B', respectively, of the membrane. A specimen is added to an extraction tube containing buffer in order to release SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens. During the test, extracted antigens, if present in the specimen, bind to anti-SARS-CoV-2 antibodies and/or anti-influenza virus type A and/or B antibodies conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-SARS-CoV-2 antibodies as well as by anti-influenza virus type A and B antibodies in the test line regions 'T' as well as 'A' and 'B', respectively. Excess coloured particles are captured in each control line region 'C'. The presence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a positive result. The absence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a negative result.

The formation of a coloured line in each control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

## 4. Reagents and Materials Supplied

- 20 NADAL® COVID-19 Ag+Influenza A/B test cassettes\*
- Additional material provided according to 93/42/EEC:  
Due to possible supply shortages of COVID-19 related accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

### a) 20 sterile swabs, CE0086

 Puritan Medical Products Company LLC  
31 School Street  
Guilford, Maine 04443-0149 USA (authorised EU representative EMERGO EUROPE, The Hague, The Netherlands)

### b) 20 sterile swabs, CE0197

 Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,  
Jiangsu 225109 China (authorised EU representative Llins Service & Consulting GmbH,  
Obere Seegasse 34/2, 69124 Heidelberg,  
Germany)

## c) 20 sterile swabs, CE0197

 CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxiu Road, Haimen, Jiangsu province  
(authorised EU representative WellKang Ltd,  
16 Castle St, Dover, CT16 1PW, UK)

## d) 20 sterile swabs, Copan Floqswabs; CE0123

 Copan Italia S.p.A.,  
Via Perotti 10, 25125 Brescia, Italy

## e) 20 sterile swabs, CE0123

 Zhejiang Gongdong Medical Technology Co., Ltd  
No.10 Beiyuan Ave., Huangyan, 318020 Taizhou,  
Zhejiang, People's Republic of China (authorised  
EU representative Shanghai International Holding  
Corp. GmbH (Europe), Eiffestraße 80, 20537  
Hamburg, Germany)

- 20 extraction tubes incl. dropper caps
- 2 buffer bottles (7 mL each)\*\*
- 1 reagent holder
- 1 package insert

\*containing the preservative sodium azide: <0.1%

\*\*Buffer containing the following preservatives: sodium azide: <0.1 mg/mL

The detergents contained in the buffer lyse and neutralise the virus.

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below exemption threshold.

**5. Additional Materials Required**

- Timer

**6. Storage & Stability**

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

**7. Warnings and Precautions**

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.

- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear appropriate personal protective equipment, such as face mask, isolation gown, gloves and eye protection during specimen collection, preparation and test procedure.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Further specimen processing and patient management should follow local guidelines and regulations.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

**8. Specimen Collection and Preparation****Nasopharyngeal specimen:**

- Insert the swab into the nostril, parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.
- Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.
- Slowly remove the swab while rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

**Note:**


- Use only synthetic fibre swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit further testing.
- Swab specimens should be tested immediately after collection. Use freshly collected specimens for best test performance.
- If not tested immediately, swab specimens can be stored in viral transport media (VTM) without denaturing agents at 2-8°C for 24 hours after collection.
- Viral transport media (VTM) without denaturing agents as well as VTM containing virus-inactivating agents and evaluated with the test can be used during specimen storage for the subsequent antigen detection using the NADAL® COVID-19 Ag+Influenza A/B Test. In order to influence the sensitivity as little as possible, a low volume of VTM (max. 1 mL). The test procedure when using viral transport media (VTM) is available separately:
  - under the following link:  
<https://www.nal-vonminden.com/en/vtm-ifu>

- by request at +49 941 29010-0 or info@nal-vonminden.com
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.

## 9. Test Procedure

**Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.**


1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.
2. Gently mix the buffer by carefully swiveling the bottle.
3. Holding the buffer bottle vertically and without touching the edge of the tube, add 10 drops to the extraction tube.
4. Insert the swab with the collected specimen into the tube. Swirl the swab and squeeze it 10-15 times by compressing the wall of the extraction tube against the swab to extract the antigens contained in the swab.
5. Let the solution stand for 2 minutes.
6. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.
7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
8. Place the test cassette on a clean and level surface.
9. Attach a dropper cap to the extraction tube, invert the tube and transfer 2 drops of the extracted solution to each specimen well (S) of the test cassette.
10. Start the timer.
11. Wait for the coloured line(s) to appear. Read the test result after exactly 15 minutes. Do not interpret the result after more than 15 minutes.

## 10. Result Interpretation

### Positive for COVID-19 Ag:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'T'.

### Positive for influenza virus type A:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'A'.

### Positive for influenza virus type B:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'B'.


**Note:** A combination of positive results for various parameters is possible.

**Note:** The colour intensity in the test line region 'T', 'A' or 'B' may vary depending on the concentration of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in the specimen. Any shade of colour in the test line region 'T', 'A' or 'B' should be considered a positive result. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.


### Negative:

A coloured line develops in each control line region 'C'. No lines develop in the test line regions 'T', 'A' and 'B'.

### Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

## 11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line developing in each control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.*

## 12. Limitations

- The NADAL® COVID-19 Ag+Influenza A/B Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens only. Neither the quantitative value nor the rate of increase/decrease in the concentration of SARS-CoV-2 as well as influenza virus type

A and B viral nucleoprotein antigens can be determined with this qualitative test.

- The NADAL® COVID-19 Ag+Influenza A/B Test only detects the presence of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in specimens and should not be used as the sole criterion for a diagnosis of COVID-19 and influenza A or B.
- Both viable and non-viable SARS-CoV-2 viruses as well as influenza viruses type A and B can be detected using the NADAL® COVID-19 Ag+Influenza A/B Test.
- The sections 'Specimen Collection and Preparation' as well as 'Test Procedure' must be followed closely while testing. Failure to follow them may lead to inaccurate test results because the antigen concentration in the swab is highly dependent on the correct procedure.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- In the course of SARS-CoV-2 or influenza A or B infection, the concentration of viral nucleoprotein antigens may fall below the detection limit of the test.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a SARS-CoV-2 or influenza A or B infection and should be confirmed via molecular assay.

### 13. Expected Values

SARS-CoV-2 viral particles are normally present in the respiratory tracts of COVID-19 patients. A positive test result can indicate an acute infection. Virus concentrations in nasopharyngeal swab specimens may vary in the course of the disease and might fall below the detection limit of rapid tests, even though patients are still showing symptoms. Conversely, the virus might continue be detectable over long periods of time even in convalescent patients. Possible infectiousness of test subjects cannot be ruled out based on negative test results.

On average, 3 to 10% of unvaccinated adults get infected with influenza each year (the number of infections varies significantly depending on the currently most prevalent strains). A positive rapid test result can indicate the presence of influenza virus type A and/or B viral nucleoprotein antigens in the specimen material, suggesting an acute infection with influenza A and/or B. Since patients might also acquire co-infections with other severe respiratory diseases, additional testing for other viral and bacterial infectious agents is recommended.

### 14. Performance Characteristics

#### Clinical performance

##### Diagnostic sensitivity and specificity

The NADAL® COVID-19 Ag+Influenza A/B Test was evaluated with clinical specimens whose status was confirmed using RT-PCR ( $C_t$  range positive: 20-37). The sensitivity for SARS-CoV-2 was calculated for the range from high to medium viral load ( $C_t$  20-30) and from high to very low viral load ( $C_t$  20-37).

The results are presented in the following tables.

RT-PCR, $C_t$ 20-30			
NADAL® COVID-19 Ag+Influenza A/B Test (COVID-19 Ag)	Positive	Negative	Total
Positive	120	0	120
Negative	3	161	164
Total	123	161	284

Diagnostic sensitivity ( $C_t$  20-30): 97.6% (93.1% - 99.2%)\*

Overall agreement ( $C_t$  20-30): 98.9% (96.9% - 99.6%)\*

Diagnostic specificity:

>99.9% (97.7% - 100%)\*

\*95% confidence interval

RT-PCR, $C_t$ 20-37			
NADAL® COVID-19 Ag+Influenza A/B Test (COVID-19 Ag)	Positive	Negative	Total
Positive	150	0	150
Negative	37	161	198
Total	187	161	348

Diagnostic sensitivity ( $C_t$  20-37): 80.2% (73.9% - 85.3%)\*

Overall agreement ( $C_t$  20-37): 89.4% (85.7% - 92.2%)\*

Diagnostic specificity:

>99.9% (97.7% - 100%)\*

\*95% confidence interval

In order to show the strong dependence of the diagnostic sensitivity on the viral load, the following table demonstrates the SARS-CoV-2 sensitivity for different  $C_t$  value ranges of the reference PCR:

$C_t$ range	Sensitivity
20 – 25	97.12%
20 – 30	97.56%
20 – 32	96.21%
20 – 35	85.71%
20 – 37	80.21%

Please note that  $C_t$  values may vary between different PCR systems at the same virus concentration.

A total of 190 nasopharyngeal swabs from patients with suspected influenza A viral infection were tested using the NADAL® COVID-19 Ag+Influenza A/B Test. Of these specimens, 78 were confirmed to be positive and 112 were confirmed to be negative by cell culture.

The results are presented in the following table:

Cell culture			
NADAL® COVID-19 Ag+Influenza A/B Test (Influenza A)	Positive	Negative	Total
Positive	65	9	74
Negative	13	103	116
Total	78	112	190

Positive agreement with cell culture: 83.3% (73.5%-90%)\*

Negative agreement with cell culture: 92% (85.4% - 95.7%)\*

Total agreement with cell culture: 88.4% (83.1% - 92.2%)\*

\*95% confidence interval

A total of 190 nasopharyngeal swabs from patients with suspected influenza B viral infection were tested using the NADAL® COVID-19 Ag+Influenza A/B Test. Of these specimens, 85 were confirmed to be positive and 105 were confirmed to be negative by cell culture.

The results are presented in the following table:

NADAL® COVID-19 Ag+Influenza A/B Test (Influenza B)	Cell culture			
		Positive	Negative	Total
	Positive	70	6	76
	Negative	15	99	114
	Total	85	105	190

Positive agreement with cell culture: 82.4% (72.9% - 89%)\*

Negative agreement with cell culture: 94.3% (88.1% - 97.4%)\*

Total agreement with cell culture: 88.9% (83.7% - 92.7%)\*

\*95% confidence interval

### Analytical performance

#### Analytical sensitivity study (detection limit)

The detection limit of the NADAL® COVID-19 Ag+Influenza A/B Test is  $2 \times 10^{2.4}$  TCID<sub>50</sub>/mL and was determined with a SARS-CoV-2 control with a known virus titre.

The detection limit of the NADAL® COVID-19 Ag+Influenza A/B Test is 0.4 ng/mL for recombinant SARS-CoV-2 nucleoprotein.

Concentrations determined as the detection limit for each influenza A/B strain tested are listed below:

Influenza A: A2/Aichi/2/68(H3N2),  $2.3 \times 10^3$  CEID<sub>50</sub>/test\*

Influenza B: Hong Kong 5/72,  $3.5 \times 10^3$  CEID<sub>50</sub>/test\*

\*CEID<sub>50</sub>: Chicken Embryo Infectious Dose

### Interfering substances

The following substances, normally present in respiratory specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below and showed no interference with the NADAL® COVID-19 Ag+Influenza A/B Test.

Substance	Concen- tra- tion	Substance	Concen- tra- tion
3 OTC nasal sprays	10%	Guaiacol glyceryl ether	20 mg/mL
3 OTC mouthwashes	10%	Mucin	1%
3 OTC sore throat liquids	10%	Mupirocin	250 µg/mL
4-acetamido-phenol	10 mg/mL	Oxymetazoline	10 mg/mL
Acetylsalicylic acid	20 mg/mL	Phenylephrine	10 mg/mL
Albuterol	20 mg/mL	Phenylpropanamine	20 mg/mL
Chlorpheniramine	5 mg/mL	Relenza® (zanamivir)	20 mg/mL
Dexamethasone	5 mg/mL	Rimantadine	500 ng/mL
Dextromethorphan	10 mg/mL	Tamiflu® (oseltamivir)	100 mg/mL
Diphenhydramine	5 mg/mL	Tobramycin	40 mg/mL
Doxylamine succinate	1 mg/mL	Triamcinolone	14 mg/mL
Flunisolide	3 mg/mL		

### Cross-reactivity

Specimens spiked with the following pathogens were tested using the NADAL® COVID-19 Ag+Influenza A/B Test (left side of the test cassette: COVID-19 Ag Test):

HCoV-HKU1, HCoV-OC43, HCoV-NL63, HCoV-229E, measles virus, *Streptococcus pneumoniae*, Epstein-Barr virus, *Bordetella parapertussis*, influenza A (H1N1) pdm09, influenza A (H3N2), influenza A (H5N1), influenza A (H7N9), influenza A (H7N7), influenza B Victoria lineage, Influenza B Yamagata lineage, *Haemophilus influenzae*, *Candida albicans*, *Mycobacterium tuberculosis*, respiratory syncytial virus, adenovirus, parainfluenza virus type 1, 2, 3, human metapneumovirus, rhinovirus, coxsackievirus type A16, norovirus, mumps virus, *Legionella pneumophila*, *Mycoplasma pneumoniae*, *Chlamydia pneumoniae*, *Streptococcus pyogenes*, *Streptococcus agalactiae*, group C *Streptococcus*, *Staphylococcus aureus*.

No cross-reactivity with the specimens was observed when tested using the NADAL® COVID-19 Ag+Influenza A/B Test (left side of the test cassette: COVID-19 Ag Test).

Specimens spiked with the following pathogens were tested using the NADAL® COVID-19 Ag+Influenza A/B Test (right side of the test cassette: Influenza A/B Test):

Human adenovirus B, C; adenovirus type 10, 18; coxsackievirus type A9, B5; human herpesvirus type 2, 5; herpes simplex virus type 1; human rhinovirus type 2, 14, 16; mumps virus; Murine respiriviruses; respiratory syncytial virus; Rubella virus; SARS-CoV-2, human coronavirus OC43; echovirus type 2, 3, 6; measles virus; parainfluenza virus type 2, 3; varicella-zoster virus.

*Acinetobacter calcoaceticus*, *Bacteroides fragilis*, *Bordetella pertussis*, *Bacillus cereus*, *Bacillus subtilis*, *Bordetella parapertussis*, *Moraxella catarrhalis*, *Chlamydia pneumoniae*, *Corynebacterium diphtheriae*, *Citrobacter freundii*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Gardnerella vaginalis*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Lactobacillus casei*, *Lactobacillus plantarum*, *Legionella pneumophila*, *Listeria monocytogenes*, *Mycobacterium avium*, *Mycobacterium intracellulare*, *Mycobacterium tuberculosis*, *Mycoplasma pneumoniae*, *Neisseria meningitidis*, *Neisseria sicca*, *Neisseria subflava*, *Nocardia asteroides*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Serratia liquefaciens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, groups A, B, C, F, G *Streptococcus*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus salivarius*, *Streptococcus sanguis*, *Yersinia enterocolitica*.

### Precision

### Repeatability and reproducibility

Precision for SARS-CoV-2 and Influenza A/B was established by testing 10 replicates of negative, low positive and high positive controls.

Reproducibility for SARS-CoV-2 was established by testing triplicates of negative, low and high positive controls. Testing was performed by 3 operators using 3 independent NADAL® COVID-19 Ag+Influenza A/B test lots at 3 different sites on 5 separate days.

Reproducibility for Influenza A/B was established by testing 10 replicates of negative, low and high positive controls. Testing was performed using 3 independent NADAL® COVID-19 Ag+Influenza A/B test lots.

The NADAL® COVID-19 Ag+Influenza A/B Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

## 15. References

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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#### 15. Referencias

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181–192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, TrendsMicrobiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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a) 20 tamponi sterili, CE0086

Puritan Medical Products Company LLC  
31 School Street  
Guilford, Maine 04443-0149 USA (rappresentante autorizzato EU EMERGO EUROPE, The Hague, The Netherlands)

b) 20 tamponi sterili, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,

Jiangsu 225109 China (rappresentante autorizzato EU Llins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)  
c) 20 tamponi sterili, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxiu Road, Haimen, Jiangsu province (rappresentante autorizzato EU WellKang Ltd, 16 Castle St, Dover, CT16 1PW, UK)

d) 20 tamponi sterili, Copan Floqswabs; CE0123  
 Copan Italia S.p.A., Via Perotti 10, 25125 Brescia, Italy

e) 20 tamponi sterili, CE0123

Zhejiang Gongdong Medical Technology Co., Ltd  
No.10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, People's Republic of China (rappresentante autorizzato EU Shanghai International Holding Corp. GmbH (Europe), Eiffestraße 80, 20537 Hamburg, Germany)

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#### 15. Bibliografia

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, TrendsMicrobiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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#### 15. Referências

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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## 15. Reference

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181–192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weis S, Leibowitz JL. Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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## 15. Lähteet

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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a) 20 sterila provtagningspinnar, CE0086

Puritan Medical Products Company LLC  
31 School Street

Guilford, Maine 04443-0149 USA (auktoriserad EU-representant EMERGO EUROPE, The Hague, The Netherlands)

b) 20 sterila provtagningspinnar, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou, Jiangsu 225109 China (auktoriserad EU-representant Lins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

## c) 20 sterila provtagningspinnar, CE0197

 CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxiu Road, Haimen, Jiangsu province  
(auktoriserad EU-representant WellKang Ltd, 16 Castle St, Dover, CT16 1PW, UK)

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## 15. Referenser

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *TrendsMicrobiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: *Fields, Virology*, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

Rev. 1, 2020-10-12 KS

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	in-vitro-Diagnostika	in-vitro diagnostic medical device	Dispositif médical de diagnostic in-vitro	Producto sanitario para diagnóstico in-vitro	Dispositivo medico-diagnóstico in-vitro	Tylko do diagnostyki in-vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkity	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käytöohjeita	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøyne
	Dispositivo médico para diagnóstico in-vitro	Diagnostický zdravotnický prostředek in-vitro	in-vitro - diagnostikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in-vitro-diagnostik	Medisch hulpmiddel voor in-vitro-diagnostiek	Medicinsk udstyr til in-vitro-diagnostik	in-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensning
	Código de lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får ikke återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotrebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogové číslo	Luetelonumero	Listnummer	Catalogus nummer	Best il i ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumääärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrekkelig for <n> tester

**Our Teams****Germany:****Regensburg**

Tel: +49 941 290 10-0  
 Fax: +49 941 290 10-50

**Moers**

Tel: +49 2841 99820-0  
 Fax: +49 2841 99820-1

**Austria:**

Tel: +49 941 290 10-29  
 Free Tel: 0800 291 565  
 Fax: +49 290 10-50  
 Free Fax: 0800 298 197

**UK & Ireland:**

Tel: +49 941 290 10-18  
 Free Tel –UK: 0808 234 1237  
 Free Tel – IRE: 1800 555 080  
 Fax: +49 290 10-50

**France:**

France Tel: 0800 915 240  
 France Fax: 0800 909 493

**Switzerland**

Swiss Tel: 0800 564 720  
 Swiss Fax: 0800 837 476

**Belgium**

Belgium Tel: 0800 718 82  
 Belgium Fax: 0800 747 07

**Luxembourg**

Lux, Tel: 800 211 16  
 Lux, Fax: 800 261 79

**Spain:**

Tel: +49 941 290 10-759  
 Free Tel: 900 938 315  
 Fax: +49 941 290 10-50  
 Free Fax: 900 984 992

**Italy:**

Tel: +49 941 290 10-34  
 Fax: +49 941 290 10-50

**Poland:**

Tel: +49 941 290 10-44  
 Free Tel: 00 800 491 15 95  
 Fax: +49 941 290 10-50  
 Free Fax: 00 800 491 15 94

**Portugal:**

Tel: +49 941 290 10-735  
 Tel, Verde: 800 849 230  
 Fax: +49 941 290 10-50  
 Fax Verde: 800 849 229

**Netherlands:**

Tel: +31 30 75 600  
 Free Tel: 0800 0222 890  
 Fax: +31 70 30 30 775  
 Free Fax: 0800 024 9519

**Nordic countries:****Denmark**

Tel: +31 703075 605  
 Free Tel: 808 887 53

**Finland**

Tel: +31 703075 606  
 Free Tel: 0800 918 263  
 Free Fax: 0800 918 262

**Norway**

Tel: +31 703075 605  
 Free Tel: 800 16 731

**Sweden**

Tel: +31 703075 605  
 Free Tel: 020 79 09 06

**Laboratory Diagnostics Team:**

Tel: +49 941 290 10-40  
 Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany  
[www.nal-vonminden.com](http://www.nal-vonminden.com) • [info@nal-vonminden.com](mailto:info@nal-vonminden.com)  
 Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH